


Contact: Annie
 Email: annie@aliance-china.com.cn
 Tel: +86 21 6252 5669
 Mobile: +86 135 8576 5009
 Skype: Annieding1986

Specification		
Image		
Name	Disposable Respirator	
Respirator Filter Class	N95	
Model	L-288	
Size	Universal	
Exhalation Valve	NO	
Model	Non-Sterile	
Color	White	
Packing	20pcs/box, 20boxes/carton	
Theoretical Techical Data		
Filtration Efficiency	BFE ≥ 95%	BFE(Bacterial filtration efficiency)
Certification	NIOSH N95 approved	
Feature & Benefits		
*Multi-material composed structure, layer by layer protection, effectively isolate dust, degrade harmful substances and prevent ubiquitous pollution in the air		
*Three dimensional protection		
*Comfortable nose clip for perfect fit protection		
*Head Strap for comfortable wearing		
* Provide at least 95% filtration efficiency against certain non-oil based particles		
*Suitable for personal daily health protection: these masks are designed to keep away smoke, debris, residue and dusts from you		
*Wide applications: can be applied in electronic manufacturing, dust-free workshops, catering services, food processing, schools, motorcycles, beauty, environmental cleaning, public occasions and many other uses		
Precautions, Warnings & Tips		
* Check the integrity of the package before use. It is strictly forbidden to use it when package damage is found.		
* The product is provided non-sterile, and it cannot be used in an environment with aseptic requirements without sterilization.		
* Not to be used with beard or other facial hair that will prevent seal		
* Place Respirator on your face so that the foam rests on your nose and the bottom panel opens under chin. Pull the top strap over your head and position it high on the back of your head. Pull the bottom strap over your head and position it around the neck and below the ears.		